Identifier:

Subfolder:

Folder:



CAPT_PRG_0100Cl1 Point of Care Urine Pregnancy (Technologist Choice) Training Checklist

CAPT_PRG_0100Cl1

POCT PREG TEST Type:

Type: Checklist (2yr) Effective on: 2025-02-05

Point of Care Urine Pregnancy (Technologist Choice) Training Checklist CAPT_PRG_200Cl1

Version #:

| 1 | TASK/COMPETENCY | COMPLETED |
|--|--|-----------|
| 1. | Complete Learning Hub: BCC-Vancouver Centre Point of Care Test - | |
| | hCG One Step Pregnancy Test (35064) course and Quiz | |
| 2. | Review the Technologist Choice One-Step Pregnancy Test Procedure. | |
| 3. | Demonstrates hand hygiene protocol and wears gloves when performing | |
| | both the quality control and patient sample. | |
| 4. | Bio-Rad Urine Quality Controls: Provided by the lab monthly. | |
| | Stored at room temperature. | |
| | Date open and expired on tubes – good for 31 days after opening at | |
| | room temperature. | |
| 5. | Perform Quality Control : | |
| | Frequency: Performed weekly or when opening new lot of kits | |
| | Enter date, lot number and expiration date of kit. | |
| | Levels 1 (blue/negative) & Level 2 (red/positive) | |
| | Open 2 kits, add 3 drops from each control to the test well. Set | |
| | timer for 5 minutes. | |
| | Document results on POCT Pregnancy QC Form: | |
| | o Pass or Fail | |
| | o Initial or Operator ID | |
| 6. | Labels Urine Sample Containers with TWO patient identifiers . | |
| 7. | Understand Sample requirement and sample stability protocols - including | |
| | sample settling time for cloudy specimens. | |
| 8. | Labels test cassette with TWO patient identifiers . | |
| 9. | Dispense three full drops urine in sample area | |
| 10. | Set timer for 5 minutes. Do not read past 10 minutes! | |
| 11. | Looks for coloured line in C and T region at 5 minutes – know test | |
| | limitations. | |
| | Positive – Two coloured lines appear | |
| | Negative – One coloured line control region © and no coloured line | |
| | in test region (T) | |
| | Invalid – No coloured line appears in the Control region © | |
| 12. | Record result in patient chart. | |
| 13. | Understand ongoing competency requirements: 1 set of QC per year and | |
| | completion of Learning Hub Course & Quiz. | |
| Employee: Signature: | | |
| | , I understand the significance of and am able to perform all the above items a | nd I will |
| follow applicable procedures. | | |
| Trained by (Lab/CNE/CNC): Date: | | |
| For further information, contact Laboratory Point of Care Laboratory | | |
| Local: 672482 email: poctlabbcca@bccancer.bc.ca | | |
| Scan and email completed checklist to: poctlabbcca@bccancer.bc.ca | | |

Written by: Reviewed by: Reviewed on: Renewed by:

Renewed on:

Ronny Garbuio

Approved by (sign.):

Approved by (name):
Approved on:
Revision Date:

Gang Wang 2025-01-31 2027-02-05

Documents used outside of OMNI are uncontrolled.